

# COUNTY OF SAN LUIS OBISPO HEALTH AGENCY PUBLIC HEALTH LABORATORY

Phone: (805) 781-5507 | Fax: (805) 781-1023 www.sloPublicHealth.org/Lab

#### **SAMPLE COLLECTION & TRANSPORT INSTRUCTIONS**

#### **Heterotrophic Plate Count Water Test**

8045

#### **PAYMENT OPTIONS**

- 1. Payment for water testing can be made by cash, check, or credit card at the reception desk upon sample drop-off. The reception desk may be closed from 12 to 1 pm on business days.
- 2. When payment is not made at sample drop-off, an invoice will be mailed to you. The invoice can be paid by cash, check, or credit card, and instructions will appear with the invoice.
- 3. For any billing questions, please contact (805) 781-5507 (option 1).
- 4. For fee information, refer to the laboratory website (www.sloPublicHealth.org/Lab).

#### **SPECIMEN STORAGE & TRANSPORT**

Temperature: Refrigerated 2-8°C (36-46°F), not frozen.

<u>Maximum Holding Time</u>: 8 hours (time from collection to test set-up in laboratory incubator). <u>Exceptions</u>: for any exceptions to collection, storage, or transport, please call lab first to discuss.

#### KIT (BOTTLE) STORAGE

<u>Temperature</u>: room temperature 15–25°C (59–77°F).

Expiration: bottle must be used prior to expiration date listed on kit label.

#### **ACCEPTABLE SPECIMEN TYPES**

Fill Line

Water: 100.0 mL in Idexx Bottle

Water level must reach the black fill line or be slightly over.

A sample under the line will be rejected.

#### **TEST ORDER OPTIONS FOR THIS KIT**

8045: Heterotrophic Plate Count (Drinking, Dental, Spa, and

Swimming Pool water)

Other: Call Laboratory to discuss

# Idexx Water Collection Bottle

#### **REJECTION CRITERIA**

- **Holding Time Exceeded:** Holding time greater than 8 hours.
- No Sample ID on Bottle or Sample ID Mismatch.
- Sample Volume Not Sufficient: Sample volume less than 100 mL.
- Inappropriate Holding Temperature: Temperature upon arrival is less than 0°C or contains ice. Note: The recommended storage temperature is refrigerated (2–8°C.). There is no maximum temperature, but the samples are to be kept chilled during transport and storage.

• Sample Leaked in Transit: Container lid not properly threaded.

Last Updated: July 3, 2025 Page 1 of 2

## SPECIMEN COLLECTION & TRANSPORT INSTRUCTIONS

### Bacteriologic (Total Coliforms/E. coli) Water Test

8040, 8350

- **Inappropriate Sample Container:** Sample submitted in container other than those provided by the lab. Note: our containers are assessed to meet the required quality control parameters.
- Evidence of Excess Chlorination: Sample has residual chlorine based on a strip test.
- **Sample Appears Overly Discolored or Contains Excess Sediment:** If the water contains a small amount of precipitate or background color, it may be tested.

#### **SAMPLE COLLECTION**

1. Use **only sterile containers supplied by this laboratory**. Sample bottles can be picked up and dropped off at any of the three Public Health Department locations below.

	Pick-Up Sample Bottles	<b>Drop-Off Samples</b>
SLO*	M-F 8-5	M-Th** by 3:30 p.m.
Paso*	M-F 8-5	M-Th** by 11:00 a.m.

<sup>\*</sup>Closed on County holidays

- 2. Do not collect samples from a swivel head faucet. When taking water samples, remove the aerator; **flush the lines** by opening the cold water tap and running water for 2 to 3 minutes.
- 3. Carefully open the bottle **do not touch the inside** of the bottle or lid.
- 4. **Do not rinse** out the bottle and **do not touch or dump** out the white powder.
- 5. Fill the bottle carefully to **the 100 ml line** on the bottle (**DO NOT UNDERFILL** or OVERFILL), making sure the powder stays in the bottle. Close the cap tightly.
- 6. **Fill in all the requested fields** on the laboratory requisition (test request form). For guidance filling out the requisition, see Environmental Testing Requisition Instructions document or contact us. Create a sample ID number and write it on both the bottle and the requisition.
- 7. Water must begin cooling, to refrigerator temperature, immediately after collection. **Use cooler** with ice packs or equivalent for transport. Do not freeze.
- 8. Cost of testing is posted at http://www.sloPublicHealth.org/Lab, payable when the water sample is delivered to any Public Health Department. If testing shows the drinking water is unsafe to drink, a "Retest" sample will incur the same fee.
- 9. Please note on the laboratory requisition if the water test is for FHA, VA, or Cal-Vet Loans. We are unable to test for metals or minerals currently.

For any assistance, call the Laboratory at (805) 781-5508 or contact Public Health in your area:

San Luis Obispo Health Dept. Paso Robles Health Dept. (805) 781-5508 (805) 237-3050
2191 Johnson Ave 805 4<sup>th</sup> Street, Suite #201
SLO, 93401 Paso Robles, 93446

Last Updated: July 3, 2025 Page 2 of 2

<sup>\*\*</sup>Water samples can be accepted on Fridays on an emergency basis, call Laboratory for approval.